


IR-FTIR ACCESSORIES

Pellet Die Set



Pellet Die Set (Circle)

- Hardened Stainless Steel
- Highly Polished pellets
- Evacuatable for sample pellet clarity and quality

Die Size	6mm	10mm	13mm	20mm	25mm	32mm	40mm	50mm
Part No.:	DS5	DS10	DS13	DS20	DS25	DS32	DS40	DS50
Part No. Total Height	65	65	65	65	95	95	95	95
Plunger Height	37	37	37	37	37	1.37	65	65
Base Diameter	8	8	8	8	13	13	13	13
Single Anvil Height	27	0.25"	27	27	27	27	27	27
Price	Nil	Nil	Nil	Nil	Nil	Nil	Nil	Nil

KBR Press

5Ton | 10Ton

KBR Press Manual



Model	HP 15T	HP 25T	HP 40T
Pressure Range	0-15T	0-25T	0-10T
Piston Diameter	56 mm	85mm	45mm
Maximum Piston Stroke	20mm	25mm	20mm
Table Diameter	95 mm	100mm	90mm
Workspace	100 X 130mm	120x150mm	90x120mm
Pressure Stability	1MPa/5min	1MPa/5min	1MPa/5min
Weight	45kg	65kg	40kg
Exterior Size	250 X 350mm	280 X 350 mm	250X350mm
Price	Nil	Nil	Nil