



Features:

Large LCD screen for easy operation in the dark, visual comfort

Touch screen operation, testing easier

High-precision insulation sensor, collect data and temperature reliable

Using halogen heating and drying, rapid heating, the test time is short

Function:

Use for the rapid determination free water of chemicals, grain, mineral, biological products, food, pharmaceutical raw materials, paper, textile raw materials and other types of samples.

Model	<i>sml-16A</i>
Max capacity	20g or 50g or 100g
readability	1mg
Moisture measurement accuracy	$\pm 0.5\%$ (sample $\geq 2g$)
Moisture content determination readability	0.02%~0.1% (sample $\geq 2g$)
Temperature control adjustment tolerance	$\pm 1^{\circ}\text{C}$
Interface	USB
Heating temperature range	60~160 $^{\circ}\text{C}$
Heating time range	0~99 min
Heating	Halogen lamp
Analytical method	Manual, timed, automatic
Moisture range	0%~100%
Dimensions	350x215x192mm
Net weight	4.5kg



Halogen moisture meter

Uniform heating/rapid heating/accurate testing

- ✓ Color touch screen
- ✓ USB interface
- ✓ Rapid drying
- ✓ Easy to operate



Shreeval
Medical Lab
Industries

*Unit 59 /B1 Building 5Th Floor/ Plot No. C -39/A
Gami Industrial Park / TTC Industrial Area
Pawane/MIDC
Navi Mumbai / 400 705 / India
9967282110*